

## Product datasheet for **TP725913**

### HBZ Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Hemoglobin subunit zeta/HBZ(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-R142
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Hemoglobin subunit zeta/HBZ is produced by E.coli expression system and the target gene encoding S2-R142 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3050
UniProt ID:	<a href="#">P02008</a>



[View online »](#)