

## Product datasheet for **TP725921**

### IL6 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Interleukin-6/IL-6(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	P29-M212
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Interleukin-6/IL-6 is produced by E.coli expression system and the target gene encoding P29-M212 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3569
UniProt ID:	<a href="#">P05231</a>



[View online »](#)