

Product datasheet for **TP725922**

SLC1A4 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Neutral amino acid transporter A/SLC1A4(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G439-L532
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Neutral amino acid transporter A/SLC1A4 is produced by E.coli expression system and the target gene encoding G439-L532 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6509
UniProt ID:	P43007
Synonyms:	Alanine/serine/cysteine/threonine transporter 1; ASCT-1; Solute carrier family 1 member 4
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)