

Product datasheet for **TP725928**

Ataxin 3 (ATXN3) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ataxin-3/ATXN3(Q291-Q299del)(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-K370
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Ataxin-3/ATXN3 is produced by E.coli expression system and the target gene encoding E2-K370 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4287
UniProt ID:	P54252
Synonyms:	Machado-Joseph disease protein 1; Spinocerebellar ataxia type 3 protein
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)