

Product datasheet for **TP725935**

Annexin A11 (ANXA11) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Annexin A11(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G344-D505
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Annexin A11 is produced by E.coli expression system and the target gene encoding G344-D505 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	311
UniProt ID:	P50995
Synonyms:	56 kDa autoantigen; Annexin XI; Annexin-11; Calyculin-associated annexin 50; CAP-50



[View online »](#)