

## Product datasheet for **TP725937**

### Malignant T cell amplified sequence 1 (MCTS1) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MCT-1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	F2-E165
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human MCT-1 is produced by E.coli expression system and the target gene encoding F2-E165 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	28985
UniProt ID:	<a href="#">Q9ULC4</a>
Synonyms:	Malignant T-cell-amplified sequence 1; Multiple copies T-cell malignancies
Protein Families:	Druggable Genome



[View online »](#)