

Product datasheet for **TP725942**

TAF11 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Transcription initiation factor TFIID subunit 11/TAF11(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	D2-F211
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Transcription initiation factor TFIID subunit 11/TAF11 is produced by E.coli expression system and the target gene encoding D2-F211 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	6882
UniProt ID:	Q15544



[View online »](#)