

## Product datasheet for **TP725944**

### Nuclear Factor 1 (NFIA) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NF1-A(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Y2-A200
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human NF1-A is produced by E.coli expression system and the target gene encoding Y2-A200 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4774
UniProt ID:	<a href="#">Q12857</a>
Synonyms:	Nuclear factor 1 A-type; Nuclear factor 1/A; CCAAT-box-binding transcription factor; CTF;Nuclear factor I/A; NF-I/A; NFI-A; TGGCA-binding protein
Protein Families:	Transcription Factors



[View online »](#)