

Product datasheet for **TP725963**

ADA Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Adenosine deaminase/ADA(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S110-L363
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Adenosine deaminase/ADA is produced by E.coli expression system and the target gene encoding S110-L363 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	100
UniProt ID:	P00813
Synonyms:	Adenosine aminohydrolase
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Primary immunodeficiency, Purine metabolism



[View online »](#)