

Product datasheet for **TP725966**

AKR1B1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Aldo-keto reductase family 1 member B1/AKR1B1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E151-F274
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Aldo-keto reductase family 1 member B1/AKR1B1 is produced by E.coli expression system and the target gene encoding E151-F274 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	231
UniProt ID:	<u>P15121</u>



[View online »](#)