

Product datasheet for **TP725973**

ATPase Inhibitory Factor 1 (ATPIF1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ATPase inhibitor, mitochondrial/ATP5IF1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G26-D106
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human ATPase inhibitor, mitochondrial/ATP5IF1 is produced by E.coli expression system and the target gene encoding G26-D106 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	93974
UniProt ID:	Q9UII2



[View online »](#)