

Product datasheet for **TP725975**

RAC3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ras-related C3 botulinum toxin substrate 3/RAC3(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	V46-T190
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Ras-related C3 botulinum toxin substrate 3/RAC3 is produced by E.coli expression system and the target gene encoding V46-T190 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5881
UniProt ID:	P60763



[View online »](#)