

Product datasheet for **TP725986**

Frequenin (NCS1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NCS-1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-V190
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human NCS-1 is produced by E.coli expression system and the target gene encoding G2-V190 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	23413
UniProt ID:	P62166
Synonyms:	Neuronal calcium sensor 1; Frequenin homolog; Frequenin-like protein; Frequenin-like ubiquitous protein
Protein Families:	Druggable Genome



[View online »](#)