

Product datasheet for **TP725991**

NHP2L1 (SNU13) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NHP2-like protein 1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-V128
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human NHP2-like protein 1 is produced by E.coli expression system and the target gene encoding T2-V128 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4809
UniProt ID:	P55769
Synonyms:	High mobility group-like nuclear protein 2 homolog 1; OTK27; SNU13 homolog; hSNU13; U4/U6.U5 small nuclear ribonucleoprotein SNU13; U4/U6.U5 tri-snRNP 15.5 kDa protein; NHP2-like protein 1, N-terminally processed
Protein Pathways:	Spliceosome



[View online »](#)