

## Product datasheet for **TP725992**

### Tumor necrosis factor alpha induced protein 8 like protein 2 (TNFAIP8L2) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TNFAIP8-like protein 2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-L184
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human TNFAIP8-like protein 2 is produced by E.coli expression system and the target gene encoding E2-L184 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	79626
UniProt ID:	<a href="#">Q6P589</a>
Synonyms:	Tumor necrosis factor alpha-induced protein 8-like protein 2; TIPE2; TNF alpha-induced protein 8-like protein 2; Inflammation factor protein 20



[View online »](#)