

Product datasheet for **TP726004**

SAA4 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Serum amyloid A-4 protein/SAA4(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E19-Y130
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Serum amyloid A-4 protein/SAA4 is produced by E.coli expression system and the target gene encoding E19-Y130 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	6291
UniProt ID:	P35542



[View online »](#)