

Product datasheet for **TP726014**

TRAX (TSNAX) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Translin-associated protein X/TSNAX
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-S290
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Translin-associated protein X/TSNAX is produced by E.coli expression system and the target gene encoding M1-S290 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7257
UniProt ID:	Q99598
Synonyms:	Translin-associated factor X



[View online »](#)