

## Product datasheet for **TP726014**

### TRAX (TSNAX) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Translin-associated protein X/TSNAX
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-S290
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Translin-associated protein X/TSNAX is produced by E.coli expression system and the target gene encoding M1-S290 is expressed.
Stability:	12 months from date of despatch
Locus ID:	7257
UniProt ID:	<a href="#">Q99598</a>



[View online »](#)