

Product datasheet for **TP726021**

ATP2B4 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PMCA4/ATP2B4
Species:	Human
Expression cDNA Clone or AA Sequence:	S172-L356
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human PMCA4/ATP2B4 is produced by E.coli expression system and the target gene encoding S172-L356 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	493
UniProt ID:	P23634
Synonyms:	Plasma membrane calcium-transporting ATPase 4; Matrix-remodeling-associated protein ; Plasma membrane calcium ATPase isoform 4; Plasma membrane calcium pump isoform 4
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Calcium signaling pathway



[View online »](#)