

Product datasheet for **TP726023**

DCUN1D1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human DCN1-like protein 1 (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-V259
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human DCN1-like protein 1 is produced by E.coli expression system and the target gene encoding M1-V259 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	54165
UniProt ID:	Q96GG9
Synonyms:	DCNL1; DCUN1 domain-containing protein 1; Defective in cullin neddylation protein 1-like protein 1;Squamous cell carcinoma-related oncogene
Protein Families:	Druggable Genome



[View online »](#)