

## Product datasheet for **TP726024**

### SCGF (CLEC11A) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human CLEC11A Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	G19-R245
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human CLEC11A Protein is produced by E.coli expression system and the target gene encoding G19-R245 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6320
UniProt ID:	<a href="#">Q9Y240</a>
Synonyms:	C-type lectin superfamily member 3; Lymphocyte secreted C-type lectin; Osteolectin; Stem cell growth factor; p47;CLECSF3; LSLCL; SCGF
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)