

Product datasheet for **TP726027**

PCGF3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Polycomb group RING finger protein 3/PCGF3
Species:	Human
Expression cDNA Clone or AA Sequence:	L2-L242
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Polycomb group RING finger protein 3/PCGF3 is produced by E.coli expression system and the target gene encoding L2-L242 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10336
UniProt ID:	Q3KNV8
Synonyms:	RING finger protein 3A; RNF3; RNF3A



[View online »](#)