

Product datasheet for **TP726029**

MCFD2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MCFD2 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	E27-Q146
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human MCFD2 Protein is produced by E.coli expression system and the target gene encoding E27-Q146 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	90411
UniProt ID:	Q8NI22
Synonyms:	Multiple coagulation factor deficiency protein 2; Neural stem cell-derived neuronal survival protein; SDNSF



[View online »](#)