

## Product datasheet for **TP726031**

### CERT1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ceramide transfer protein/hCERT
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-F598
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Ceramide transfer protein/hCERT is produced by E.coli expression system and the target gene encoding S2-F598 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10087
UniProt ID:	<a href="#">Q9Y5P4</a>
Synonyms:	Collagen type IV alpha-3-binding protein; Goodpasture antigen-binding protein; GPBP; START domain-containing protein 11; StAR-related lipid transfer protein 11; CERT1; CERT; COL4A3BP; STARD11



[View online »](#)