

Product datasheet for **TP726035**

HNRPH1 (HNRNPH1) Human Recombinant Protein

Product data:

| | |
|--|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human hnRNP H |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | N220-A449 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human hnRNP H is produced by E.coli expression system and the target gene encoding N220-A449 is expressed with a His tag at the C-terminus. |
| Storage: | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 3187 |
| UniProt ID: | P31943 |
| Synonyms: | Heterogeneous nuclear ribonucleoprotein H;HNRNPH1; HNRPH1 |



[View online »](#)