

Product datasheet for **TP726036**

NKRF Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NFkB-repressing factor
Species:	Human
Expression cDNA Clone or AA Sequence:	N345-N690
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human NFkB-repressing factor is produced by E.coli expression system and the target gene encoding N345-N690 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55922
UniProt ID:	O15226
Synonyms:	NF-kappa-B-repressing factor; NRF;Protein ITBA4; NKRF; ITBA4
Protein Families:	Transcription Factors



[View online »](#)