

## Product datasheet for **TP726042**

### TRIM26 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Tripartite motif-containing protein 26/TRIM26 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	N151-W270
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Tripartite motif-containing protein 26/TRIM26 Protein is produced by E.coli expression system and the target gene encoding N151-W270 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7726
UniProt ID:	<a href="#">Q12899</a>
Synonyms:	Acid finger protein (AFP); RING finger protein 95; Zinc finger protein 173
Protein Families:	Druggable Genome



[View online »](#)