

Product datasheet for **TP726043**

GORASP2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GRS2
Species:	Human
Expression cDNA Clone or AA Sequence:	L152-P452
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human GRS2 is produced by E.coli expression system and the target gene encoding L152-P452 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	26003
UniProt ID:	Q9H8Y8
Synonyms:	Golgi reassembly-stacking protein 2; Golgi phosphoprotein 6 (GOLPH6); Golgi reassembly-stacking protein of 55 kDa (GRASP55); p59; GORASP2
Protein Families:	Druggable Genome



[View online »](#)