

Product datasheet for **TP726044**

REEP4 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Receptor expression-enhancing protein 4/REEP4 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	R157-S257
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Receptor expression-enhancing protein 4/REEP4 Protein is produced by E.coli expression system and the target gene encoding R157-S257 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	80346
UniProt ID:	Q9H6H4
Synonyms:	C8orf20; PP432
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)