

Product datasheet for **TP726046**

TNIP1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TNFAIP3-interacting protein 1
Species:	Human
Expression cDNA Clone or AA Sequence:	A287-Q636
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human TNFAIP3-interacting protein 1 is produced by E.coli expression system and the target gene encoding A287-Q636 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10318
UniProt ID:	Q15025
Synonyms:	A20-binding inhibitor of NF-kappa-B activation 1 (ABIN-1); HIV-1 Nef-interacting protein; Nef-associated factor 1 (Naf1); Nip40-1; Virion-associated nuclear shuttling protein (VAN; hVAN); TNIP1



[View online »](#)