

## Product datasheet for **TP726047**

### HSPA2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Heat shock 70 kDa protein 2
Species:	Human
Expression cDNA Clone or AA Sequence:	L290-D639
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Heat shock 70 kDa protein 2 is produced by E.coli expression system and the target gene encoding L290-D639 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3306
UniProt ID:	<a href="#">P54652</a>



[View online »](#)