

Product datasheet for **TP726048**

FAM62B (ESYT2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human E-Syt2
Species:	Human
Expression cDNA Clone or AA Sequence:	N559-L773
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human E-Syt2 is produced by E.coli expression system and the target gene encoding N559-L773 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	57488
UniProt ID:	A0FGR8
Synonyms:	Extended synaptotagmin-2; Chr2Syt; FAM62B; KIAA1228
Protein Families:	Transmembrane



[View online »](#)