

Product datasheet for **TP726053**

Tryptophanyl tRNA synthetase (WARS) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Tryptophan--tRNA ligase, cytoplasmic/WARS1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	P2-V171
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Tryptophan--tRNA ligase, cytoplasmic/WARS1 Protein is produced by E.coli expression system and the target gene encoding P2-V171 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	7453
UniProt ID:	P23381



[View online »](#)