

## Product datasheet for **TP726055**

### VAT1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Synaptic vesicle membrane protein VAT-1 homolog
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-T148
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Synaptic vesicle membrane protein VAT-1 homolog is produced by E.coli expression system and the target gene encoding S2-T148 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10493
UniProt ID:	<a href="#">Q99536</a>
Synonyms:	VAT1 Protein
Protein Families:	Druggable Genome



[View online »](#)