

Product datasheet for **TP726056**

PCBP2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Poly(rC)-binding protein 2/PCBP2 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	D2-E168
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Poly(rC)-binding protein 2/PCBP2 Protein is produced by E.coli expression system and the target gene encoding D2-E168 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5094
UniProt ID:	Q15366



[View online »](#)