

## Product datasheet for **TP726058**

### DUT Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human dUTPase
Species:	Human
Expression cDNA Clone or AA Sequence:	P2-A130
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human dUTPase is produced by E.coli expression system and the target gene encoding P2-A130 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	1854
UniProt ID:	<a href="#">P33316</a>
Synonyms:	Deoxyuridine 5'-triphosphate nucleotidohydrolase; mitochondrial; dUTP pyrophosphatase
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Pyrimidine metabolism



[View online »](#)