

Product datasheet for **TP726060**

DYNLT1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Dynein light chain Tctex-type 1/Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-S112
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Dynein light chain Tctex-type 1/Protein is produced by E.coli expression system and the target gene encoding E2-S112 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6993
UniProt ID:	P63172
Synonyms:	Protein CW-1; T-complex testis-specific protein 1 homolog



[View online »](#)