

Product datasheet for **TP726063**

CHCHD10 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human CHCHD10 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	A17-P142
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human CHCHD10 Protein is produced by E.coli expression system and the target gene encoding A17-P142 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	400916
UniProt ID:	Q8WYQ3
Synonyms:	Coiled-coil-helix-coiled-coil-helix domain-containing protein 10; mitochondrial; Protein N27C7-4



[View online »](#)