

Product datasheet for **TP726065**

PLPBP Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PLP homeostasis protein
Species:	Human
Expression cDNA Clone or AA Sequence:	N21-H275
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human PLP homeostasis protein is produced by E.coli expression system and the target gene encoding N21-H275 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	11212
UniProt ID:	O94903
Synonyms:	Pyridoxal phosphate homeostasis protein; Proline synthase co-transcribed bacterial homolog protein; Pyridoxal phosphate-binding protein



[View online »](#)