

Product datasheet for **TP726067**

LAP3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cytosol aminopeptidase/LAP3 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	I40-V280
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Cytosol aminopeptidase/LAP3 Protein is produced by E.coli expression system and the target gene encoding I40-V280 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	51056
UniProt ID:	P28838



[View online »](#)