

Product datasheet for **TP726080**

Elongin A (TCEB3) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Elongin-A/EloA
Species:	Human
Expression cDNA Clone or AA Sequence:	A679-R798
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Elongin-A/EloA is produced by E.coli expression system and the target gene encoding A679-R798 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6924
UniProt ID:	Q14241
Synonyms:	Elongin 110 kDa subunit; RNA polymerase II transcription factor SIII subunit A1; SIII p110; Transcription elongation factor B polypeptide 3
Protein Families:	Transcription Factors



[View online »](#)