

Product datasheet for **TP726081**

GDNF Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GDNF
Species:	Human
Expression cDNA Clone or AA Sequence:	S78-I211
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human GDNF is produced by E.coli expression system and the target gene encoding S78-I211 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2668
UniProt ID:	P39905
Synonyms:	Glial cell line-derived neurotrophic factor; Astrocyte-derived trophic factor (ATF)
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)