

## Product datasheet for **TP726082**

### DDX58 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Antiviral innate immune response receptor RIG-I
Species:	Human
Expression cDNA Clone or AA Sequence:	D792-K925
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Antiviral innate immune response receptor RIG-I is produced by E.coli expression system and the target gene encoding D792-K925 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	23586
UniProt ID:	<a href="#">O95786</a>



[View online »](#)