

Product datasheet for **TP726085**

FGF5 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human FGF-5
Species:	Human
Expression cDNA Clone or AA Sequence:	S145-G268(K238N,P245S)
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human FGF-5 is produced by E.coli expression system and the target gene encoding S145-G268(K238N,P245S) is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	2250
UniProt ID:	P12034



[View online »](#)