

## Product datasheet for **TP726085**

### FGF5 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human FGF-5
Species:	Human
Expression cDNA Clone or AA Sequence:	S145-G268(K238N,P245S)
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human FGF-5 is produced by E.coli expression system and the target gene encoding S145-G268(K238N,P245S) is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2250
UniProt ID:	<a href="#">P12034</a>
Synonyms:	Fibroblast growth factor 5; Heparin-binding growth factor 5 (HBGF-5); Smag-82
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	MAPK signaling pathway, Melanoma, Pathways in cancer, Regulation of actin cytoskeleton



[View online »](#)