

Product datasheet for **TP726087**

RILP Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Rab-interacting lysosomal protein
Species:	Human
Expression cDNA Clone or AA Sequence:	A73-K259
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Rab-interacting lysosomal protein is produced by E.coli expression system and the target gene encoding A73-K259 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	83547
UniProt ID:	Q96NA2
Synonyms:	RILP Protein



[View online »](#)