

Product datasheet for **TP726088**

Caspase 8 (CASP8) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human CASP-8/Caspase-8
Species:	Human
Expression cDNA Clone or AA Sequence:	S141-D496
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human CASP-8/Caspase-8 is produced by E.coli expression system and the target gene encoding S141-D496 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	841
UniProt ID:	Q14790



[View online »](#)