

## Product datasheet for **TP726092**

### PP2A-beta (PPP2CB) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PP2A-beta
Species:	Human
Expression cDNA Clone or AA Sequence:	D2-K144
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human PP2A-beta is produced by E.coli expression system and the target gene encoding D2-K144 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5516
UniProt ID:	<a href="#">P62714</a>
Synonyms:	Serine/threonine-protein phosphatase 2A catalytic subunit beta isoform
Protein Families:	Druggable Genome, Phosphatase
Protein Pathways:	Long-term depression, Oocyte meiosis, TGF-beta signaling pathway, Tight junction, Wnt signaling pathway



[View online »](#)