

Product datasheet for **TP726096**

POLE2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human DNA polymerase epsilon subunit 2
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-P200
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human DNA polymerase epsilon subunit 2 is produced by E.coli expression system and the target gene encoding A2-P200 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5427
UniProt ID:	P56282



[View online »](#)