

Product datasheet for **TP726097**

GALE Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human UDP-glucose 4-epimerase
Species:	Human
Expression cDNA Clone or AA Sequence:	G198-N340
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human UDP-glucose 4-epimerase is produced by E.coli expression system and the target gene encoding G198-N340 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	2582
UniProt ID:	Q14376



[View online »](#)